




**AUTOSTOP
FILTER**



**PRIMING
FILTER**



TRO-SOLUSET
AUTOSTOP

INFUSION SET

ADVANTAGES

of infusion set with **AUTOSTOP**



Our latest innovative feature adds extra safety and eliminates the risk of air embolism in patients receiving infusions by preventing air from passing through the **AUTOSTOP** filter.

- ▶ Specially engineered **AUTOSTOP** filter (15 micron) prevents the entry of air in the line when the IV bottle is empty.
- ▶ **AUTOSTOP** filter maintains a constant fluid level. No air can get through to the patient, thus eliminating the risk of air embolism.

Features

- ▶ The special **Priming Filter** ensures that the line fills automatically (auto prime) thus maintains a closed system until connected.

- ▶ Helps in preventing **Healthcare Associated Infections (HAI)**.



- ▶ **No more messy spillage of IV fluid** on hands, floor and bedding, also reduces wastage of IV fluid.



Secure Disposal of Spike:
Spike docking facility for safe and secure disposal.

Roller Body
Specially designed roller body with smooth operation and precise flow control, also tubing can be parked after priming for stand by procedures.



- ▶ PVC free flexible drip chamber for easy and fast priming (Drop rate 20 drops/ml)
- ▶ DEHP free Infusion line with new generation non-phthalate plasticizers, reduces phthalate related toxicity and ensures patient safety
- ▶ Universally suitable for Gravity Infusion as well as Pressure Infusions up to 2 bars

CARING FOR GLOBAL HEALTH