



TRO-SEPTO SWAB

DISINFECTION CAP

TRO-SEPTO SWAB

DISINFECTION CAP

Advantages: Safety and Comfort



TRO-SEPTO SWAB, disinfection cap, is intended for use on needle-free valves as a disinfecting cleaner prior to line access and act as a cover between line accesses. The device will disinfect the needle-free valves one minute after application and protects from contamination between accesses for up to seven days if not removed. It shall not be used as a replacement for a Combi-Stopper cap.

- Ergonomic design for smooth usage
- Provides instant active and efficient disinfection for 7 days if not removed
- Compatible with all needle free luer access systems
- Eliminates inter-personal variation in compliance of scrub the hub method of disinfection and saves nursing time
- PVC Free, Latex Free, DEHP Free
- contains 70% Isopropyl alcohol
- Gamma Sterile
- Available in strips of 10 caps.
- Single unit packs upon request

This product is marketed and distributed by TROGE MEDICAL GmbH and CE-certified by a German authorized Notified Body.



CARING FOR GLOBAL HEALTH

TROGE MEDICAL GMBH | Milchstrasse 19 | 20148 Hamburg | Germany | Tel. +49 (0) 40 / 44 18 44 0 | info@trogemedical.de

www.trogemedical.de