

TRO-VEBSITE safety **B protect**

BLOOD CONTROL IV CANNULA

IV Cannula with **B protect**

Advantages: Safety and Comfort



Our latest innovative feature:
A Silicon Septum to prevent
blood spillage



Features


- Blood control IV cannula reduces the risk of blood exposure and contamination
- Helping to maintain a clean access site
- Passive safety mechanism automatically encapsulates sharp needle bevel, preventing accidental needle stick injuries
- Needle safety guard significantly reduces the risk of transmission of blood borne pathogens like HBV, HCV & HIV
- Advanced polyurethane catheter material ensures better bio-compatibility & longer indwelling time
- Perforated wings provide ventilation & prevent skin maceration
- Reduces the risk of air-in-line when IV bottle runs dry
- Compatible with power injectors upto 300 PSI
- Latex free, PVC free, DEHP free




User friendly, self activating safety mechanism prevents needle stick injuries



Cannula reduces the risk of blood exposure and contamination



Flexible perforated soft wings



Also available: Quick rapid flash needle technology ensures successful venipuncture

Also available with injection port:



Needle Gauge	Catheter O.D. (mm)	Catheter Length (mm)	Flow Rate (ml/min)	Flow Rate (ml/sec)	Ref. No. (non-ported)	Ref. No. (ported)
18 G ■	1.3	45	65	10	90822	90834
18 G ■	1.3	32	70	10	90823	90835
20 G ■	1.1	32	47	7	90824	90836
22 G ■	0.9	25	28	6	90825	90837
24 G ■	0.7	19	17	4	90826	90838